

## ALTRONIX MODELO TPX-III pH-Bench Top

# Neocientífica

## ALTRONIX®



### Model: TPX-III

#### Feature:

1. Microprocessor based designed, and rugged housing with splash proof keyboard.
2. The large backlight LCD display pH and Temperature, or millivolts simultaneously.
3. Automatic or manual temperature compensation and degree °C/°F switchable.
4. All reading feature stamp meeting GLP guidelines.
5. Stores up to 64 reading which you can download data to your computer via RS232 output.

### Specifications :

	TPX-III		
	pH	mV	Temp.
Range:	0-14.00	±2000	-5-105°C
Accuracy:	±0.01 + 1 digit	0.1% mV/F.S.	±0.2°C + 1 digit
Resolution:	0.01	1	0.1°C
Compensation:	ATC 0-100°C		
Communication:	RS 232		
Storage data	64 reading		
Power:	DC 9V by AC adaptor		
Dimensions:	200 x 160 x 65 mm		
Weight:	600 g (complete kit)		

### Accessories for complete kit:

1. pH electrode(U.S.A.).
2. Temperature probe(Pt 1000).
3. AC/DC adaptor.
4. Electrode Holder.
5. RS232 software & connection cable.
6. Instruction manual.